



Product Datasheet

Anti-CLDN6 Reference Antibody (orb1805484)


Description	Anti-CLDN6 Reference Antibody) is expressed from CHO. The heavy chain type is IgG1, and the light ch...
Reactivity	Human
Endotoxins	< 1 EU/mg
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, Kinetics
Target	CLDN6
Preservatives	100 mM Pro-Ac 20 mM Arg pH 5.0
Storage	-80°C
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	> 95%
MW	150 kDa
Expiration Date	12 months from date of receipt.



Anti-CLDN6
Reference
Antibody
(DS-
9606a)...



The purity
of Anti-
CLDN6
Reference
Antib...



Immobilized
human
CLDN6 VLP
Protein at
2...