


Product Datasheet

Anti-AA2AR / Adenosine A2aR Reference Antibody (orb1805493)


Description

Anti-AA2AR / Adenosine A2aR Reference


	chain type is IgG1, ...
Reactivity	Human
Endotoxins	< 1 EU/mg
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, Kinetics
Target	AA2AR / Adenosine A2aR
Preservatives	100 mM Pro-Ac 20 mM Arg pH 5.0
Storage	-80°C
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	> 95%
MW	146.66 kDa
Expiration Date	12 months from date of receipt.



Anti-AA2AR
/ Adenosine
A2aR
Reference
An...



The purity
of Anti-
AA2AR /
Adenosine
A2a...



Immobilized
human
A2aR VLP
at 2 µg/mL
C...