

Product Datasheet

Anti-TrkB / NTRK2 Reference Antibody (orb1805521)

Description Anti-TrkB / NTRK2 Reference Antibody is expressed from CHO. The heavy chain type is hulgG2, and the ...

Reactivity Human

Endotoxins < 1 EU/mg

Conjugation Unconjugated

Tested Applications ELISA, FA, FACS, Kinetics

Target TrkB / NTRK2

Preservatives 100 mM Pro-Ac 20 mM Arg pH 5.0

Storage -80°C

Note For research use only


Isotype IgG

Clonality Monoclonal


Purity > 95%

MW 145.5 kDa

Expiration Date 12 months from date of receipt.



Anti-TrkB /
NTRK2
Reference
Antibody
(Ri...



The purity
of Anti-TrkB
/ NTRK2
Referenc...