

## Product Datasheet

### Anti-GAD65 Reference Antibody (orb1805537)

**Description** Anti-GAD65 Reference Antibody is expressed from CHO. The heavy chain type is hulgG1, and the light C...

**Reactivity** Human

**Endotoxins** < 1 EU/mg

**Conjugation** Unconjugated

**Tested Applications** ELISA, FA, FACS, Kinetics

**Target** GAD65

**Preservatives** 100 mM Pro-Ac 20 mM Arg pH 5.0

**Storage** -80°C

**Note** For research use only


**Isotype** IgG

**Clonality** Monoclonal


**Purity** > 95%

**MW** 145.5 kDa

**Expiration Date** 12 months from date of receipt.



Anti-GAD65  
Reference  
Antibody  
(U.Washing...



The purity of  
Anti-GAD65  
Reference  
Antib...