

## Product Datasheet

### Anti-Siglec-2 / CD22 Reference Antibody (orb1805578)


|                            |                                                                                                         |
|----------------------------|---------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Anti-Siglec-2 / CD22 Reference Antibody is expressed from CHO. The heavy chain type is hulgG1, and t... |
| <b>Reactivity</b>          | Human                                                                                                   |
| <b>Endotoxins</b>          | < 1 EU/mg                                                                                               |
| <b>Conjugation</b>         | Unconjugated                                                                                            |
| <b>Tested Applications</b> | ELISA, FA, FACS, Kinetics                                                                               |
| <b>Target</b>              | Siglec-2 / CD22                                                                                         |
| <b>Preservatives</b>       | 100 mM Pro-Ac 20 mM Arg pH 5.0                                                                          |
| <b>Storage</b>             | -80°C                                                                                                   |
| <b>Note</b>                | For research use only                                                                                   |
| <b>Isotype</b>             | IgG                                                                                                     |
| <b>Clonality</b>           | Monoclonal                                                                                              |
| <b>Purity</b>              | > 95%                                                                                                   |
| <b>MW</b>                  | 145.5 kDa                                                                                               |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                         |



Anti-Siglec-2 / CD22 Reference Antibody ...



The purity of Anti-Siglec-2 / CD22 Refer...



Immobilized human CD22 His at 2 µg/mL c...