

Product Datasheet

Anti-Siglec-15 / CD33L3 Reference Antibody (orb1805582)

Description Anti-Siglec-15 / CD33L3 Reference Antibody is expressed from CHO. The heavy chain type is hulgG1, an...

Reactivity Human

Endotoxins < 1 EU/mg

Conjugation Unconjugated

Tested Applications ELISA, FA, FACS, Kinetics

Target Siglec-15 / CD33L3

Preservatives 100 mM Pro-Ac 20 mM Arg pH 5.0

Storage -80°C

Note For research use only


Isotype IgG

Clonality Monoclonal

Purity > 95%

MW 145.5 kDa

Expiration Date 12 months from date of receipt.



Anti-Siglec-
15 /
CD33L3
Reference
Antibo...



The purity
of Anti-
Siglec-15 /
CD33L3
Re...