

Product Datasheet

Anti-BST2 / CD317 Reference Antibody (orb1805584)

Description Anti-BST2 / CD317 Reference Antibody is expressed from CHO. The heavy chain type is hulgG1, and the ...

Reactivity Human

Endotoxins < 1 EU/mg

Conjugation Unconjugated

Tested Applications ELISA, FA, FACS, Kinetics

Target BST2 / CD317

Preservatives 100 mM Pro-Ac 20 mM Arg pH 5.0

Storage -80°C

Note For research use only


Isotype IgG

Clonality Monoclonal


Purity > 95%

MW 145.5 kDa

Expiration Date 12 months from date of receipt.



Anti-BST2 /
CD317
Reference
Antibody
(SB...



The purity
of Anti-
BST2 /
CD317
Referenc...