

Product Datasheet

Anti-CYR61 / CCN1 Reference Antibody (orb1805593)

Description Anti-CYR61 / CCN1 Reference Antibody is expressed from CHO. The heavy chain type is hulgG1, and the ...

Reactivity Human

Endotoxins < 1 EU/mg

Conjugation Unconjugated

Tested Applications ELISA, FA, FACS, Kinetics

Target CYR61 / CCN1

Preservatives 100 mM Pro-Ac 20 mM Arg pH 5.0

Storage -80°C

Note For research use only


Isotype IgG

Clonality Monoclonal


Purity > 95%

MW 145.5 kDa

Expiration Date 12 months from date of receipt.



Anti-CYR61
/ CCN1
Reference
Antibody
(Ro...



The purity
of Anti-
CYR61 /
CCN1
Referenc...