

Product Datasheet

Anti-PDCD1 / PD-1 / CD279 Reference Antibody (orb1805601)

Description Anti-PDCD1 / PD-1 / CD279 Reference Antibody is expressed from CHO. The heavy chain type is hulgG4SP...

Reactivity Human

Endotoxins < 1 EU/mg

Conjugation Unconjugated

Tested Applications ELISA, FA, FACS, Kinetics

Target PDCD1 / PD-1 / CD279

Preservatives 100 mM Pro-Ac 20 mM Arg pH 5.0

Storage -80°C

Note For research use only


Isotype IgG

Clonality Monoclonal

Purity > 95%

MW 145.5 kDa

Expiration Date 12 months from date of receipt.



Anti-PDCD1
/ PD-1 /
CD279
Reference
Anti...



The purity
of Anti-
PDCD1 /
PD-1 /
CD279 ...