

Product Datasheet

Anti-OSMR Reference Antibody (orb1805635)

Description Anti-OSMR Reference Antibody is expressed from CHO. The heavy chain type is hulgG4SP, and the light ...

Reactivity Human

Endotoxins < 1 EU/mg

Conjugation Unconjugated

Tested Applications ELISA, FA, FACS, Kinetics

Target OSMR

Preservatives 100 mM Pro-Ac 20 mM Arg pH 5.0

Storage -80°C

Note For research use only


Isotype IgG

Clonality Monoclonal


Purity > 95%

MW 145.5 kDa


Expiration Date 12 months from date of receipt.



Anti-OSMR
Reference
Antibody
(vixarelima...



The purity of
Anti-OSMR
Reference
Antibo...



Immobilized
human
OSMR His at
2 µg/mL c...