

Product Datasheet

Anti-HGF / SF Reference Antibody (orb1805741)

Description Anti-HGF / SF Reference Antibody is expressed from CHO. The heavy chain type is hulgG1, and the ligh...

Reactivity Human

Endotoxins < 1 EU/mg

Conjugation Unconjugated

Tested Applications ELISA, FA, FACS, Kinetics

Target HGF / SF

Preservatives 100 mM Pro-Ac 20 mM Arg pH 5.0

Storage -80°C

Note For research use only


Isotype IgG

Clonality Monoclonal


Purity > 95%

MW 145.5 kDa

Expiration Date 12 months from date of receipt.



Anti-HGF /
SF
Reference
Antibody
(TAK-70...



The purity
of Anti-HGF
/ SF
Reference
An...