

Product Datasheet

Anti-TRP1 / TYRP1 Reference Antibody (orb1805826)

Description Anti-TRP1 / TYRP1 Reference Antibody is expressed from CHO. The heavy chain type is hulgG1, and the ...

Reactivity Human

Endotoxins < 1 EU/mg

Conjugation Unconjugated

Tested Applications ELISA, FA, FACS, Kinetics

Target TRP1 / TYRP1

Preservatives 100 mM Pro-Ac 20 mM Arg pH 5.0

Storage -80°C

Note For research use only


Isotype IgG

Clonality Monoclonal


Purity > 95%

MW 145.5 kDa

Expiration Date 12 months from date of receipt.



Anti-TRP1 /
TYRP1
Reference
Antibody
(fl...



The purity
of Anti-
TRP1 /
TYRP1
Referenc...