

Product Datasheet

Anti-AXL / UFO Reference Antibody (orb1805842)

Description Anti-AXL / UFO Reference Antibody is expressed from CHO. The heavy chain type is hulgG1, and the lig...

Reactivity Human

Endotoxins < 1 EU/mg

Conjugation Unconjugated

Tested Applications ELISA, FA, FACS, Kinetics

Target AXL / UFO

Preservatives 100 mM Pro-Ac 20 mM Arg pH 5.0

Storage -80°C

Note For research use only


Isotype IgG

Clonality Monoclonal


Purity > 95%

MW 145.42 kDa


Expiration Date 12 months from date of receipt.



Anti-AXL /
UFO
Reference
Antibody
(enapo...



The purity
of Anti-AXL /
UFO
Reference
A...



Immobilized
human Axl
FC at 2
µg/mL can...