


Product Datasheet

Mouse IGFBP2 ELISA Kit (orb180590)

Description	Mouse IGFBP2 ELISA Kit
Reactivity	Mouse
Range	125 pg/ml - 8,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: E35-Q305
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	32847 MW
Uniprot ID	P47877
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Mouse
IGFBP2
PicoKine
ELISA Kit
standard...