


Product Datasheet

Human IL1R1/II 1 Ri ELISA Kit (orb180591)

Description	Human IL1R1/II 1 Ri ELISA Kit
Reactivity	Human
Range	125 pg/ml - 8,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: D21-T332
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants and serum.
MW	65402 MW
Uniprot ID	P14778
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum is less than the lowest standard, 12.5 pg/ml.
Expiration Date	Please enquire.



Human
IL1R1
PicoKine
ELISA Kit
standard ...