


Product Datasheet

Anti-AXL / UFO Reference Antibody (orb1805936)


Description	Anti-AXL / UFO Reference Antibody is expressed from CHO. The heavy chain type is hulgG1, and the lig...
Reactivity	Human
Endotoxins	< 1 EU/mg
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, Kinetics
Target	AXL / UFO
Preservatives	100 mM Pro-Ac 20 mM Arg pH 5.0
Storage	-80°C
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	> 95%
MW	146.16 kDa
Expiration Date	12 months from date of receipt.



Anti-AXL /
UFO
Reference
Antibody
(tilve...



The purity
of Anti-AXL /
UFO
Reference
A...



Immobilized
human Axl
FC at 2
µg/mL can...