


Product Datasheet

Mouse MIP-1Alpha/CCL3 ELISA Kit (orb180594)

Description	Mouse MIP-1Alpha/CCL3 ELISA Kit
Reactivity	Mouse
Range	7.8 pg/ml - 500 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: A24-A92
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate).
MW	10345 MW
Uniprot ID	P10855
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is around 50 pg/ml.
Expiration Date	Please enquire.



Mouse
CCL3/MIP1
alpha
PicoKine
ELISA Kit...