


Product Datasheet

Anti-TNFRSF10B / TRAILR2 / CD262 Reference Antibody (orb1805962)


Description

Anti-TNFRSF10B / TRAILR2 / CD262 Reference


	chain type is h...
Reactivity	Human
Endotoxins	< 1 EU/mg
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, Kinetics
Target	TNFRSF10B / TRAILR2 / CD262
Preservatives	100 mM Pro-Ac 20 mM Arg pH 5.0
Storage	-80°C
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	> 95%
MW	145.5 kDa
Expiration Date	12 months from date of receipt.



Anti-TNFRSF10B / TRAILR2 / CD262 Referen...



The purity of Anti-TNFRSF10B / TRAILR2 / ...



Immobilized Cynomolgus TRAIL R2/DR5/ TNF...