

Product Datasheet

Anti-TREM1 / CD354 Reference Antibody (orb1805967)

Description

Anti-TREM1 / CD354 Reference Antibody is

Reactivity

hulG1, and the...

Human

Endotoxins

< 1 EU/mg

Conjugation

Unconjugated

Tested Applications

ELISA, FA, FACS, Kinetics

Target

TREM1 / CD354

Preservatives

100 mM Pro-Ac 20 mM Arg pH 5.0

Storage

-80°C

Note

For research use only

Isotype

IgG

Clonality

Monoclonal

Purity

> 95%


MW

145.5 kDa

Expiration Date

12 months from date of receipt.


MDs NR M R



Anti-TREM1
/ CD354
Reference
Antibody
(P...



The purity
of Anti-
TREM1 /
CD354
Referen...



Immobilized
human
TREM1 His
at 2 µg/mL

...