


Product Datasheet

Human N-Cadherin-2 CDH2 CD325 ELISA Kit (orb180597)

Description	Human N-Cadherin-2 CDH2 CD325 ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: D160-A724
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	99809 MW
Uniprot ID	P19022
Sensitivity	20 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is around 1500 pg/ml.
Expiration Date	Please enquire.



Human
Cadherin-
2/N-
Cadherin
PicoKine
ELI...