


Product Datasheet

Human ESM1/Endocan ELISA Kit (orb180598)

Description	Human ESM1/Endocan ELISA Kit
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: W20-R184
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	20095 MW
Uniprot ID	Q9NQ30
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is 200-800 pg/ml.
Expiration Date	Please enquire.



Human
ESM1/Endocan
PicoKine ELISA
Kit st...