

Product Datasheet

Anti-SCFR / c-Kit / CD117 Reference Antibody (orb1805985)

Description Anti-SCFR / c-Kit / CD117 Reference Antibody is expressed from CHO. The heavy chain type is hulgG1, ...

Reactivity Human

Endotoxins < 1 EU/mg

Conjugation Unconjugated

Tested Applications ELISA, FA, FACS, Kinetics

Target SCFR / c-Kit / CD117

Preservatives 100 mM Pro-Ac 20 mM Arg pH 5.0

Storage -80°C

Note For research use only

Isotype IgG

Clonality Monoclonal


Purity > 95%

MW 145.5 kDa

Expiration Date 12 months from date of receipt.



Anti-SCFR /
c-Kit /
CD117
Reference
Anti...



The purity
of Anti-
SCFR / c-Kit
/ CD117 ...