


Product Datasheet

Anti-GPC1 / Glypican-1 Reference Antibody (orb1805997)

Description	Anti-GPC1 / Glypican-1 Reference Antibody is expressed from CHO. The heavy chain type is mulgG1, and...
Reactivity	Human
Endotoxins	< 1 EU/mg
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, Kinetics
Target	GPC1 / Glypican-1
Preservatives	100 mM Pro-Ac 20 mM Arg pH 5.0
Storage	-80°C
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	> 95%
MW	145.5 kDa
Expiration Date	12 months from date of receipt.



Anti-GPC1 /
Glypican-1
Reference
Antibod...



The purity
of Anti-
GPC1 /
Glypican-1
Ref...