



Product Datasheet

Anti-TMEFF1 / Tomoregulin-1 Reference Antibody (orb1806002)

Description	Anti-TMEFF1 / Tomoregulin-1 Reference Antibody is expressed from CHO. The heavy chain type is hulg1...
Reactivity	Human
Endotoxins	< 1 EU/mg
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, Kinetics
Target	TMEFF1 / Tomoregulin-1
Preservatives	100 mM Pro-Ac 20 mM Arg pH 5.0
Storage	-80°C
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	> 95%
MW	145.5 kDa
Expiration Date	12 months from date of receipt.



Anti-TMEFF1
/
Tomoregulin-
1 Reference
An...



The purity of
Anti-TMEFF1
/
Tomoregulin-
...