


Product Datasheet


Anti-DLK1 Reference Antibody (orb1806008)

Description	Anti-DLK1 Reference Antibody is expressed from CHO. The heavy chain type is hulgG1, and the light ch...
Reactivity	Human
Endotoxins	< 1 EU/mg
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, Kinetics
Target	DLK1
Preservatives	100 mM Pro-Ac 20 mM Arg pH 5.0
Storage	-80°C
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	> 95%
MW	145.5 kDa
Expiration Date	12 months from date of receipt.




Anti-DLK1
Reference
Antibody
(LIV-1205)

...



The purity of
Anti-DLK1
Reference
Antibo...



Immobilized
human
DLK1.1 His
at 2 µg/mL...