


Product Datasheet

Human Nidogen-1 / Entactin / NID-1 / NID1 ELISA Kit (orb180602)

Description	Human Nidogen-1 / Entactin / NID-1 / NID1 ELISA Kit
Reactivity	Human
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: L29-K1114
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants and serum.
MW	136377 MW
Uniprot ID	P14543
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human
Nidogen-
1/Entactin/NID-
1 PicoKine ...