


Product Datasheet

Anti-TIE2 / CD202b Reference Antibody (orb1806027)

Description	Anti-TIE2 / CD202b Reference Antibody is expressed from CHO. The heavy chain type is hulgG4SP, and t...
Reactivity	Human
Endotoxins	< 1 EU/mg
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, Kinetics
Target	TIE2 / CD202b
Preservatives	100 mM Pro-Ac 20 mM Arg pH 5.0
Storage	-80°C
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	> 95%
MW	145.5 kDa
Expiration Date	12 months from date of receipt.



Anti-TIE2 /
CD202b
Reference
Antibody
(R...



The purity
of Anti-TIE2
/ CD202b
Referen...