


Product Datasheet

Anti-GPR73 / PROKR1 Reference Antibody (orb1806031)

Description	Anti-GPR73 / PROKR1 Reference Antibody is expressed from CHO. The heavy chain type is hulgG4SP, and ...
Reactivity	Human
Endotoxins	< 1 EU/mg
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, Kinetics
Target	GPR73 / PROKR1
Preservatives	100 mM Pro-Ac 20 mM Arg pH 5.0
Storage	-80°C
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	> 95%
MW	145.5 kDa
Expiration Date	12 months from date of receipt.



Anti-GPR73
/ PROKR1
Reference
Antibody
(...



The purity
of Anti-
GPR73 /
PROKR1
Refere...