


Product Datasheet

Mouse TIMP-2 ELISA Kit (orb180604)

Description	Mouse TIMP-2 ELISA Kit
Reactivity	Mouse
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: C27-P220
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	24328 MW
Uniprot ID	P25785
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:200, concentration in serum and plasma is 60-170 ng/ml.
Expiration Date	Please enquire.



Mouse TIMP-2 PicoKine ELISA Kit standard...