

Product Datasheet

Anti-ILDR2 Reference Antibody (orb1806057)

Description Anti-ILDR2 Reference Antibody is expressed from CHO. The heavy chain type is hulgG2, and the light C...

Reactivity Human

Endotoxins < 1 EU/mg

Conjugation Unconjugated

Tested Applications ELISA, FA, FACS, Kinetics

Target ILDR2

Preservatives 100 mM Pro-Ac 20 mM Arg pH 5.0

Storage -80°C

Note For research use only


Isotype IgG

Clonality Monoclonal

Purity > 95%

MW 145.5 kDa

Expiration Date 12 months from date of receipt.



Anti-ILDR2
Reference
Antibody
(bapotulim...



The purity of
Anti-ILDR2
Reference
Antib...