

Product Datasheet

Anti-CSF2Rb / CD131 Reference Antibody (orb1806062)

Description

Anti-CSF2Rb / CD131 Reference Antibody is

Reactivity

Human

Endotoxins

< 1 EU/mg

Conjugation

Unconjugated

Tested Applications

ELISA, FA, FACS, Kinetics

Target

CSF2Rb / CD131

Preservatives

100 mM Pro-Ac 20 mM Arg pH 5.0

Storage

-80°C

Note

For research use only

Isotype

IgG

Clonality

Monoclonal

Purity

> 95%


MW

145.5 kDa


Expiration Date

12 months from date of receipt.


kDa NR M R



Anti-CSF2Rb / CD131 Reference Antibody (...)



The purity of Anti-CSF2Rb / CD131 Refere...



Immobilized human CSF2RB His at 2 µg/mL...