



Product Datasheet

Anti-INHBA / Activin A Reference Antibody (orb1806065)

Description	Anti-INHBA / Activin A Reference Antibody is expressed from CHO. The heavy chain type is hulgG4SP, a...
Reactivity	Human
Endotoxins	< 1 EU/mg
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, Kinetics
Target	INHBA / Activin A
Preservatives	100 mM Pro-Ac 20 mM Arg pH 5.0
Storage	-80°C
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	> 95%
MW	145.5 kDa
Expiration Date	12 months from date of receipt.



Anti-INHBA
/ Activin A
Reference
Antibod...



The purity
of Anti-
INHBA /
Activin A
Ref...