



Product Datasheet

Anti-CB1 / CNR1 Reference Antibody (orb1806081)

Description	Anti-CB1 / CNR1 Reference Antibody is expressed from CHO. The heavy chain type is hulgG1, and the li...
Reactivity	Human
Endotoxins	< 1 EU/mg
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, Kinetics
Target	CB1 / CNR1
Preservatives	100 mM Pro-Ac 20 mM Arg pH 5.0
Storage	-80°C
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	> 95%
MW	145.5 kDa
Expiration Date	12 months from date of receipt.




Anti-CB1 /
CNR1
Reference
Antibody
(GFB-...



The purity
of Anti-CB1
/ CNR1
Reference

...



Human CB1
EGFP CHOS
cells were
stained w...