

Product Datasheet

Anti-TNFRSF13C / BAFFR / CD268 Reference Antibody (orb1806082)

Description

Anti-TNFRSF13C / BAFFR / CD268 Reference

chain type is hul...

Reactivity

Human

Endotoxins

< 1 EU/mg

Conjugation

Unconjugated

Tested

Applications

ELISA, FA, FACS, Kinetics

Target

TNFRSF13C / BAFFR / CD268

Preservatives

100 mM Pro-Ac 20 mM Arg pH 5.0

Storage

-80°C

Note

For research use only

Isotype

IgG

Clonality

Monoclonal

Purity

> 95%

MW

145.5 kDa

Expiration Date


12 months from date of receipt.



Anti-TNFRSF13C / BAFFR / CD268 Reference...



The purity of Anti-TNFRSF13C / BAFFR / C...



Immobilized human GPA33 / A33; Fc tag at...