

Product Datasheet

Anti-TNFRSF12A / TWEAKR / CD266 Reference Antibody (orb1806095)

Description Anti-TNFRSF12A / TWEAKR / CD266 Reference Antibody is expressed from CHO. The heavy chain type is hu...

Reactivity Human

Endotoxins < 1 EU/mg

Conjugation Unconjugated

Tested Applications ELISA, FA, FACS, Kinetics

Target TNFRSF12A / TWEAKR / CD266

Preservatives 100 mM Pro-Ac 20 mM Arg pH 5.0

Storage -80°C

Note For research use only


Isotype IgG

Clonality Monoclonal

Purity > 95%

MW 145.5 kDa

Expiration Date 12 months from date of receipt.



Anti-TNFRSF12A / TWEAKR / CD266 Reference...



The purity of Anti-TNFRSF12A / TWEAKR /

...