



Product Datasheet

Anti-STAB1 Reference Antibody (orb1806107)


Description	Anti-STAB1 Reference Antibody is expressed from CHO. The heavy chain type is hulgG4SP, and the light...
Reactivity	Human
Endotoxins	< 1 EU/mg
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, Kinetics
Target	STAB1
Preservatives	100 mM Pro-Ac 20 mM Arg pH 5.0
Storage	-80°C
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	> 95%
MW	145.5 kDa
Expiration Date	12 months from date of receipt.



Anti-STAB1
Reference
Antibody
(bexmarili...



The purity
of Anti-
STAB1
Reference
Antib...



Immobilized
human
STAB His at
2 µg/mL c...