

Product Datasheet


Anti-SIRPa / CD172a Reference Antibody (orb1806109)

Description


Anti-SIRPa / CD172a Reference Antibody is

	hulG4SP, and ...
Reactivity	Human
Endotoxins	< 1 EU/mg
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, Kinetics
Target	SIRPa / CD172a
Preservatives	100 mM Pro-Ac 20 mM Arg pH 5.0
Storage	-80°C
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	> 95%
MW	145.5 kDa
Expiration Date	12 months from date of receipt.


kDa N M R



Anti-SIRPa /
CD172a
Reference
Antibody
(...



The purity
of Anti-
SIRPa /
CD172a
Refere...



Immobilized
human
SIRPα V1
FC at4
µg/m...