



Product Datasheet

Anti-SEMA4D / CD100 Reference Antibody (orb1806110)


Description	Anti-SEMA4D / CD100 Reference Antibody is expressed from CHO. The heavy chain type is huIgG4SP, and ...
Reactivity	Human
Endotoxins	< 1 EU/mg
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, Kinetics
Target	SEMA4D / CD100
Preservatives	100 mM Pro-Ac 20 mM Arg pH 5.0
Storage	-80°C
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	> 95%
MW	145.5 kDa
Expiration Date	12 months from date of receipt.



Anti-SEMA4D /
CD100
Reference
Antibody (...)



The purity of
Anti-SEMA4D /
CD100 Refere...



Immobilized
human
SEMA4D/CD100;
His Tag ...