



Product Datasheet

Anti-PVRIG Reference Antibody (orb1806118)


Description	Anti-PVRIG Reference Antibody is expressed from CHO. It has a predicted MW of 145.5 kDa.
Reactivity	Human
Endotoxins	< 1 EU/mg
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, Kinetics
Target	PVRIG
Preservatives	100 mM Pro-Ac 20 mM Arg pH 5.0
Storage	-80°C
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	> 95%
MW	145.5 kDa
Expiration Date	12 months from date of receipt.



Anti-PVRIG
Reference
Antibody
(COM701)
O...



The purity of
Anti-PVRIG
Reference
Antib...



Immobilized
human
PVRIG His at
2 µg/mL ...