


Product Datasheet

Rat N-Cadherin-2 CDH2 CD325 ELISA Kit (orb180613)

Description	Rat N-Cadherin-2 CDH2 CD325 ELISA Kit
Reactivity	Rat
Range	0.78 ng/ml - 50 ng/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: D160-A724
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	99686 MW
Uniprot ID	Q9Z1Y3
Sensitivity	20 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is less than the lowest standard, 15.6 pg/ml.
Expiration Date	Please enquire.



Rat
Cadherin-
2/N-
Cadherin
PicoKine
ELISA...