

Product Datasheet

Anti-MASP2 Reference Antibody (orb1806148)

Description Anti-MASP2 Reference Antibody is expressed from CHO. The heavy chain type is hulgG4SP, and the light...

Reactivity Human

Endotoxins < 1 EU/mg

Conjugation Unconjugated

Tested Applications ELISA, FA, FACS, Kinetics

Target MASP2

Preservatives 100 mM Pro-Ac 20 mM Arg pH 5.0

Storage -80°C

Note For research use only


Isotype IgG

Clonality Monoclonal


Purity > 95%

MW 145.5 kDa


Expiration Date 12 months from date of receipt.



Anti-MASP2
Reference
Antibody
(narsoplim...



The purity
of Anti-
MASP2
Reference
Antib...



Immobilized
human
MASP 2A His
at 2 µg/m...