


Product Datasheet

Mouse JAM-A / F11R / Junctional Adhesion Molecule 1 ELISA Kit (orb180616)

Description	Mouse JAM-A / F11R / Junctional Adhesion Molecule 1 ELISA Kit
Reactivity	Mouse
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: K27-G238
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	32424 MW
Uniprot ID	O88792
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is around 250 pg/ml.
Expiration Date	Please enquire.



Mouse
JAM-
A/F11R
PicoKine
ELISA Kit
stan...