


Product Datasheet

Human TFF1/PS2 ELISA Kit (orb180618)

Description	Human TFF1/PS2 ELISA Kit
Reactivity	Human
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: E25-F84
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	9150 MW
Uniprot ID	P04155
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is around 50 pg/ml.
Expiration Date	Please enquire.



Human TFF1
PicoKine
ELISA Kit
standard c...