

Product Datasheet

Anti-HBEGF Reference Antibody (orb1806193)

Description Anti-HBEGF Reference Antibody is expressed from CHO. It has a predicted MW of 145.5 kDa.

Reactivity Human

Endotoxins < 1 EU/mg

Conjugation Unconjugated

Tested Applications ELISA, FA, FACS, Kinetics

Target HBEGF

Preservatives 100 mM Pro-Ac 20 mM Arg pH 5.0

Storage -80°C

Note For research use only


Isotype IgG

Clonality Monoclonal

Purity > 95%

MW 145.5 kDa

Expiration Date 12 months from date of receipt.



Anti-HBEGF
Reference
Antibody
(U3-1565) ...



The purity
of Anti-
HBEGF
Reference
Antib...