


Product Datasheet

Human TFF3/ITF ELISA Kit (orb180620)

Description	Human TFF3/ITF ELISA Kit
Reactivity	Human
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: E22-Y80
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) , saliva and urine.
MW	8641 MW
Uniprot ID	Q07654
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is 1-2 ng/ml.
Expiration Date	Please enquire.



Human TFF3
PicoKine
ELISA Kit
standard c...