


Product Datasheet

Human Midkine ELISA Kit (orb180621)

Description	Human Midkine ELISA Kit
Reactivity	Human
Range	78 pg/ml - 5,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: K23-D143
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	15585 MW
Uniprot ID	P21741
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human
Midkine
PicoKine
ELISA Kit
standar...