


Product Datasheet

Human FOLR1/Alpha Fr ELISA Kit (orb180623)

Description	Human FOLR1/Alpha Fr ELISA Kit
Reactivity	Human
Range	46.9 pg/ml - 3,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: R25-S234
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) , saliva, urine and human milk.
MW	29819 MW
Uniprot ID	P15328
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is 341-1240 pg/ml.
Expiration Date	Please enquire.



Human
FOLR1
PicoKine
ELISA Kit
standard ...