

Product Datasheet

Anti-CLEC12A / CD371 Reference Antibody (orb1806237)

Description Anti-CLEC12A / CD371 Reference Antibody is expressed from CHO. The heavy chain type is hulgG1, and t...

Reactivity Human

Endotoxins < 1 EU/mg

Conjugation Unconjugated

Tested Applications ELISA, FA, FACS, Kinetics

Target CLEC12A / CD371

Preservatives 100 mM Pro-Ac 20 mM Arg pH 5.0

Storage -80°C

Note For research use only


Isotype IgG

Clonality Monoclonal


Purity > 95%

MW 145.5 kDa

Expiration Date 12 months from date of receipt.



Anti-
CLEC12A /
CD371
Reference
Antibody ...



The purity
of Anti-
CLEC12A /
CD371
Refer...