


Product Datasheet


Anti-CLDN6 Reference Antibody (orb1806238)

Description	Anti-CLDN6 Reference Antibody is expressed from CHO. It has a predicted MW of 145.5 kDa.
Reactivity	Human
Endotoxins	< 1 EU/mg
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, Kinetics
Target	CLDN6
Preservatives	100 mM Pro-Ac 20 mM Arg pH 5.0
Storage	-80°C
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	> 95%
MW	145.5 kDa
Expiration Date	12 months from date of receipt.




Anti-CLDN6
Reference
Antibody
(IMAB027)

...



The purity
of Anti-
CLDN6
Reference
Antib...



Immobilized
human
CLDN6 VLP
at 8 µg/mL

...