

Product Datasheet

Anti-CD79b Reference Antibody (orb1806247)

Description Anti-CD79b Reference Antibody is expressed from CHO. The heavy chain type is hulgG1, and the light C...

Reactivity Human

Endotoxins < 1 EU/mg

Conjugation Unconjugated

Tested Applications ELISA, FA, FACS, Kinetics

Target CD79b

Preservatives 100 mM Pro-Ac 20 mM Arg pH 5.0

Storage -80°C

Note For research use only


Isotype IgG

Clonality Monoclonal


Purity > 95%

MW 145.5 kDa


Expiration Date 12 months from date of receipt.



Anti-CD79b
Reference
Antibody
(polatuzum...)



The purity of
Anti-CD79b
Reference
Antib...



Immobilized
human
CD79b His
at 16 µg/mL

...